

EUROTUBO® TUBES WITH VALVE CAP FOR HUMAN VENOUS BLOOD EXTRACTION

Made of transparent polypropylene, shock proof. Synthetic rubber caps, with a valve that allows the introduction and extraction of the blood sample without the need of decapping. Labelled tubes with indication of code, batch number, expiry date, which allows product traceability. Once the blood is extracted, it is introduced into the corresponding tubes, depending on the analysis to be performed **WHOLE BLOOD, SERUM or PLASMA**.

Tubes for **SERUM** obtention:

These tubes contain an additive which accelerates the clotting with gel.

Instructions for use:

- Remove the needle and discard.
- Introduce the blood sample into the tube across the valve cap.
- Invert gently from 4 to 6 times with the hand.
- Let it rest for 30 minutes vertically.
- Centrifugate at 1500 xg for 15 minutes. The accelerator eases the coagulation and the gel forms a barrier between the clot and the serum.

Tubes for **PLASMA** obtention:

Depending on the analysis to be performed the following additives are available: Lithium Heparin, EDTA dipotassium and tripotassium and Sodium citrate.

The citrates (VSG and coagulation) are in liquid form, according with the established standards.

Instructions for use:

- Remove the needle and discard.
- Introduce the blood sample into the tube across the valve cap to the label mark.
- Invert gently from 4 to 6 times with the hand.
- For total blood analysis, the test can be done immediately. For plasma test, the sample must rest or be centrifugated gently. Please see General Catalogue where more detailed and specific information is available.

Store tubes far from light in a well-ventilated area and away from any heat or ignition point.

Avoid extreme temperatures.

Do not use if the product has changed its original colour.











Do not use if dirtiness is visible.


Do not use if the expiry date is exceeded.

Do not use if broken or with missing components.

Disposal as biological waste.

Symbol glossary:

 Use-by-date	 Do not re-use	 Quantity
 Batch code	 Catalogue number	 Unique device identifier
 In Vitro Diagnostic Medical Device	 Consult instructions for use on the website www.deltalab.es/eifus	 CE marking
		 Manufacturer

 In case of a serious incident* related to the product, notify to Deltalab, S.L. as well as the competent authority of the State in which the user is established.
 *A "serious incident" is understood as one that entails the death, or serious deterioration of the health of the patient or user or a serious threat to public health.