



en

## BED PAN AND URINE COLLECTION BOTTLE

**For the collection of human biological fluids (urine and feces). Clean and disinfect before use. No sterilised.**  
**The product is not manufactured with The product is not made of natural rubber latex or dried natural rubber.**  
**Medical device.**











### Intended use:

The bed pan is intended for the collection of urine and/or feces in bedridden patients or with reduced mobility. The urine collection bottle is intended for urine collection in bedridden male patients. Both the pan and the bottle are Intended for hospital and/or home use.

### General precautions:

1. Not suitable for any other application than its intended use.
2. The use of the pan and bottle presents no difficulty for the intended end users.
3. Do not use the pan or bottle if the package is damaged.
4. It should be cleaned with detergent and disinfected with bleach or other disinfectants, before each use. The pan and bottle are suitable for use in disinfectant washers and for autoclaving.
5. The life of the pan is determined by wear and deterioration due to use. Discard when colour changes, roughness, deformations, etc. are detected.
6. Store away from heat sources.

### Symbol glossary:

 <b>REF</b>	Catalogue number	 <b>LOT</b>	Batch code	 <b>CE</b>	CE marking
 <b>MD</b>	Medical Device	 <b>QTY</b>	Quantity		Use-by-date
	Keep away from sunlight		Manufacturer		Do not use if Package is damaged
	Consult instructions for use on the website <a href="http://www.deltalab.es">www.deltalab.es</a> or the link <a href="http://bit.ly/eifus">http://bit.ly/eifus</a>				

In case of a serious incident\* related to the product, notify to Deltalab, S.L. as well as the competent authority of the State in which the user is established. \*A "serious incident" is understood as one that entails the death, or serious deterioration of the health of the patient or user or a serious threat to public health.



REVISION DATE: 2022/01



**DELTALAB, S.L.**

Plaza de La Verneda 1, Pol. Ind. La Llana,  
08191 Rubí, Barcelona. España/Spain.  
[deltalab@deltalab.es](mailto:deltalab@deltalab.es) - [www.deltalab.es](http://www.deltalab.es)

TEC 1597/0